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Storage:
The Product should be stored between 0° C- 22° C.

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Instructions for use

Patient target group
All persons for whom individual trays are required for treatment.

Intended user
Dentist

Composition Qualitative data:
Organic urethane acrylate, vinyl ester resin, 2,2-Dimethoxy-1,2-diphenylethane-1-one or by 2,2-Dimethoxy-2-phenylacetone, (+)-2,3-Dioxo-1,7,7-trimethyl-bicyclo(2,2,1)heptan-2,3-dione, Diethyl(2,4,6-trimethylbenzoyl)phosphinoxid, 2-Hydroxy-2-methyl-1-phenylpropan-1-one, 2-Ethylhexyl-4-(dimethylamino)benzoate, 2,2-Diethoxyacetophenone, soda-lime glass, amorphous silica, dyes (Pigment Red 146 = C12485 / Pigment Blue 29 = C17007 / D&C Red 30 = C17360 / Pigment Blue C17416) / Pigment Blue C17416 / Pigment Violet 19 = C17390), Aromastoffe (Optamint).

Quantitative data:
In 100 g there are 19 g acrylates/methacrylates and 80 g glass.

Intended purpose of the product
smardent Löffelmaterial is a light-curing tray material in form of prefabricated, baseplates to make rigid, custom impression trays and bite-blocks, as well as other indications. smardent Löffelmaterial can be cured in all ultraviolet and Halogen light curing units. The baseplates will be stored in a reusable light protection box that can be used for transporting or storing models, etc.

Indications
Custom impression tray

Contraindications
The material contains methacrylates. Do not use in patients with an allergy to methacrylates or other substances contained in the product.

Performance characteristics of the product
smardent Löffelmaterial polymerises with light to a solid and compact polymer.

Specific precautions
Avoid any unshielded eyes and skin and the uncured material, i.e. also with the brush. Because of the dust development during the grinding, always use a dust extractor when trimming the tray!

Preparation/Manufacture
Mark the outline of the future tray on the model.

Unterstützende Bereiche am Modell ausblenden.
To create a spacer for the impression material, block out the entire model (e.g. with a sheet of paper) and cut out stops.

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Zeiten:
Working at daylight 20 minutes

Polymerisation in ultraviolet or Halogen light-curing unit, depending on the unit's use:

1. Underside of the tray, off the model	3-5 minutes
2. Underside of the tray, off the model	3-5 minutes

Trimming:
Trim the custom tray with tungsten carbide burs and smooth with sandpaper.

Hygine:
Immerse impressions taken with smart Löffelmaterial Löffeln in a disinfectant solution based on glutaraldehyde. The time depends on the disinfectant manufacturer's directions.

The finished medical product is used by the dentist, e.g. for individual trays with impression material.

Vorsicht:
The product has been developed for use in the dental sector and must be used according to the directions. The manufacturer accepts no liability for damage that results from use other than the above. Furthermore, the user is obliged under his own responsibility to examine the material for its suitability and its usage possibilities for the intended applications, especially if the latter are not detailed in the directions for use.

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Alle serious incidents related to the device shall be reported to the manufacturer and to the competent authority of the Member State where the user and/or the patient is established.

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Composition Qualitative data:
Organic urethane acrylate, vinyl ester resin, 2,2-Dimethoxy-1,2-diphenylethane-1-one or by 2,2-Dimethoxy-2-phenylacetone, (+)-2,3-Dioxo-1,7,7-trimethyl-bicyclo(2,2,1)heptan-2,3-dione, Diethyl(2,4,6-trimethylbenzoyl)phosphinoxid, 2-Hydroxy-2-methyl-1-phenylpropan-1-one, 2-Ethylhexyl-4-(dimethylamino)benzoate, 2,2-Diethoxyacetophenone, verre sodocalcique, silice amorphe, colorants (pigment rouge 146 = C12485 / pigment bleu 29 = C17007 / D&C rouge 30 = C17360 / pigment bleu C17416) / pigment violet 19 = C17390), aromastoffes (Optamint).

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Composizione Informativa Qualitativa:
Organic urethane acrylate, vinyl ester resin, 2,2-Dimethoxy-1,2-diphenylethane-1-one or by 2,2-Dimethoxy-2-phenylacetone, (+)-2,3-Dioxo-1,7,7-trimethyl-bicyclo(2,2,1)heptan-2,3-dione, Diethyl(2,4,6-trimethylbenzoyl)phosphinoxid, 2-Hydroxy-2-methyl-1-phenylpropan-1-one, 2-Ethylhexyl-4-(dimethylamino)benzoate, 2,2-Diethoxyacetophenone, vetro sodocalcico, silice amorfa, coloranti (Pigmento Rosso 146 = C12485 / Pigmento Blu 29 = C17007 / D&C Rosso 30 = C17360 / Pigmento Blu C17416) / Pigmento Violetto 19 = C17390), aromastoffes (Optamint).

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Samenstelling kwalitatieve gegevens:
Organic urethane acrylate, vinyl ester resin, 2,2-Dimethoxy-1,2-diphenylethane-1-one or by 2,2-Dimethoxy-2-phenylacetone, (+)-2,3-Dioxo-1,7,7-trimethyl-bicyclo(2,2,1)heptan-2,3-dione, Diethyl(2,4,6-trimethylbenzoyl)phosphinoxid, 2-Hydroxy-2-methyl-1-phenylpropan-1-one, 2-Ethylhexyl-4-(dimethylamino)benzoate, 2,2-Diethoxyacetophenone, soda-lime glass, amorphous silica, dyes (Pigment Red 146 = C12485 / Pigment Blue 29 = C17007 / D&C Red 30 = C17360 / Pigment Blue C17416) / Pigment Violet 19 = C17390), aromastoffes (Optamint).