

Technical File TF15




Technical specification

Device group: Endodontic files

Doc-id: TF37.10_SPEC

Revision: 01

Pgs: 1 (2)

Technical specification											
Device trade name(s)	Sendoline® K-Reamer This specification applies to any other brand of model KR.										
Device model	KR										
Illustration of file type	 Note: The image is prior to mounting of the stopper on the file.										
File type illustration	 Reamer type K.										
Identification symbol	 Reamer type K.										
Variants/dimensions	<table border="1"><thead><tr><th><u>Lengths</u></th><th><u>Sizes</u></th></tr></thead><tbody><tr><td>21 mm</td><td>06-140</td></tr><tr><td>25 mm</td><td>06-140</td></tr><tr><td>28 mm</td><td>06-140</td></tr><tr><td>31 mm</td><td>06-140</td></tr></tbody></table> The colour of the handle is standardised as follows: 06 – pink 08 – grey 10 – purple 15 – white 20 – yellow 25 – red 30 – blue 35 – green 40 – black 45 – white 50 – yellow 55 – red 60 – blue 70 – green 80 – black 90 – white 100 – yellow	<u>Lengths</u>	<u>Sizes</u>	21 mm	06-140	25 mm	06-140	28 mm	06-140	31 mm	06-140
<u>Lengths</u>	<u>Sizes</u>										
21 mm	06-140										
25 mm	06-140										
28 mm	06-140										
31 mm	06-140										

Technical File TF15


Technical specification

Device group: Endodontic files



Doc-id: TF37.10_SPEC

Revision: 01

Pgs: 2 (2)

	110 – red 120 – blue 130 – green 140 – black
Packages	Individual sizes in blister: 6 pcs per blister of same length and size. Assortments in blister: 6 pc per blister of same length and 6 different sizes: <u>21 mm:</u> 1 x 15, 1 x 20, 1 x 25, 1 x 30, 1 x 35, 1 x 40 1 x 45, 1 x 50, 1 x 55, 1 x 60, 1 x 70, 1 x 80 <u>25 mm:</u> 1 x 15, 1 x 20, 1 x 25, 1 x 30, 1 x 35, 1 x 40 1 x 45, 1 x 50, 1 x 55, 1 x 60, 1 x 70, 1 x 80 <u>28 mm:</u> 1 x 15, 1 x 20, 1 x 25, 1 x 30, 1 x 35, 1 x 40 1 x 45, 1 x 50, 1 x 55, 1 x 60, 1 x 70, 1 x 80 <u>31 mm:</u> 1 x 15, 1 x 20, 1 x 25, 1 x 30, 1 x 35, 1 x 40 1 x 45, 1 x 50, 1 x 55, 1 x 60, 1 x 70, 1 x 80
Package illustration	 <p>Note: The image is only intended to demonstrate package design. It may include files of another model.</p>
Product description and features	<ul style="list-style-type: none">- Type 1: standard instrument (taper = 0.2) according to EN ISO 3630-1.- Shaping and cleaning instrument. File type KR: open spiral design with three cutting edges.- 60° tip.- Cutting ability. For canal enlargement.- Handheld: should be used with a screwing movement.- ISO sizes: colour coded handle.- Size and type reamer K symbol printed on handle.- Excellent debris transportation.- Handle with grooves.- Silicone stopper (white) mounted on the instrument.- Non-sterile package.- Operative part of stainless steel.- Silicone handle.

Revision history

Revision	Date	Prepared by	Reviewed by	Approved by	Remark
01	2023-01-24	Carina Hjorth, RA	Elisabet Jakobsson, Product Manager 	Henric Karsk, CEO 	Initial document.